
## **GENTLE FIX F.C CLIP**



GENTLE FIX F.C CLIP Foley Catheter Fixation Clip Solution for babies, elderly, and patients with sensitive skin.



### Why use GENTLE FIX F.C CLIP

A medical device for fixating urinary catheters at the catheter bifurcation.

#### **Features and Benefits**

- A safe use on skin
- Suitable for those whom with weakened immunity such as children and adults, and patients with serious illness because the film was made of medical silicone adhesive which minimized skin irritation.
- It can be removed without itching, pain, and damaging the skin because its adhesion lasts from former to latter part.
- It reduces the risk of trauma and pain caused by the catheter pulling or repeatedly rubbing on the bladder neck and urethra.

#### Excellent fixation

The medical silicone adhesion of film was maximized by applying COBUILD's unique technology to the film so it can be fixed stably even on sweaty skin for long time without irritation.

#### Waterproof

The film serves as a bacterial and viral barrier protecting the wound from potential contamination and provides waterproof qualities by applying polyurethane film to it so it can stay on watery and sweaty skin.

#### Natural movement

It applies revolvable clip so catheter moves in natural unison with the body.

#### COBUILD, Inc.


COBUILD B/D, 49-40, Haengju-ro 15beon-gil, Deogyang-gu, Goyang-si, Gyeonggi-do, Korea

TEL: +82-31-971-0010 FAX: +82-31-971-1112

# **GENTLE FIX F.C CLIP**



COBUILD Inc. research, manufacture and market world's finest silicone based medical products



### GENTLE FIX F.C CLIP Ordering Information

| <b>Product Code</b> | Size         | Pieces/Box |
|---------------------|--------------|------------|
| 0901                | 11.5cm x 9cm | 10/1       |